

# PATRIOT



**Item# F-12PPL-HD-AT1**  
**Item# F-12PPL-HD-AT2**

The Patriot Lift is the only completely portable aquatic lift that is independently verified to meet ADA guidelines. When deployed and ready for use, the Patriot Lift features a specifically designed seat allowing for easy lateral transfers from most wheelchairs. This gives full freedom of access to the individual. The Patriot maintains a 375 pound weight capacity.



U.S. Patent Nos.  
7,249,386 B2 &  
D507,769S

## ACCENT COLORS

Available in a variety of powder coat colors



Available in a variety of accent colors



## FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 375 lb (170 kg) weight capacity
- Setback range for Fixing Kit 26" - 31" (66 - 79 cm)
- Fixing Kit included
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty


## OPTIONS/ACCESSORIES:

F-120PPC	Cover, standard lockable, blue
F-120PPC-G	Cover, standard lockable, gray
F-120PPC-T	Cover, standard lockable, tan
F-019CA	Cycle attachment
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
F-043SCH	Solar charging station


**PROUDLY SOLD BY:**

# THE PATRIOT LIFT

## PRODUCT SPECIFICATIONS



## APPLICABLE POOL PROFILES



To confirm specific application go to [www.aquacreek.com](http://www.aquacreek.com) for our Pool Lift Selector  
Become a member of our new Dealer and Architect page. Access to CAD blocks,  
CSI specs and sales materials.