

Pride® Lifts

Commander 450

The Commander 450 exterior powered boom style lift is the perfect user-friendly lift option for pickup trucks that provides fast loading and unloading of your mobility device. Standard with a 450 lb (204.1 kg) lifting capacity for heavier power wheelchairs and scooters, it is available in 2-axis and 3-axis. This lift delivers maximum durability, with the ability to lift and lower your mobility device to the ground over the truck's tailgate.

FEATURES

- 450 lb (204.1 kg) lifting capacity
- 2-axis design: powered up/down and powered rotation left/right
- 3-axis design: powered up/down, powered rotation left/right, powered adjustable telescoping boom
- Standard single seat post or 4-point seat post carabiner kit
- Coastal community powder coated
- Telescoping mounting brackets
- Complete weather-proof design
- Waterproof push-button hand control
- Compatible with Power Wheelchairs and Scooters

BASE WEIGHT ¹	109 lb (49.4 kg)
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	70" (177.8 cm)
MOTOR	12-volt DC
WARRANTY	3-Year Transferable
VEHICLE PLATFORM	Full-Size Vans and Trucks

1. May not include some accessories
2. May require vehicle modification
3. Hitch weight carrying capacity varies by vehicle. Call for details.

