

CarryLite Transit Chair 200H

Designed for easy storage, the Evac+Chair® CarryLite Transit Chair folds neatly in half and can be stored in the smallest of areas. The chair unfolds in one swift movement for immediate use. to spiral stairs and where space is restricted.



Features

1 Lightweight aluminum construction

Increased strength and robustness for continued use

2 Integrated handles (no pinch points)

For effective grip when lifting and transporting the passenger

3 Large wheels

For increased stability when transporting over rough terrain

4 Unique locking system

When deployed it automatically locks chair into it's useable position

5 Adjustable patient straps

For increased passenger safety

6 Footrest

For passenger comfort during transportation



Includes

- LIFETIME WARRANTY
- Passenger Safety Belt

Benefits

- FDA Registered
- Easy to use
- Unfolds in one swift movement
- Compact storage
- Lightweight
- Reduces potential injury and claims
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act



Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	23"	9"	11"	16lbs	440lbs

VIDEO LINK: [Evac-Chair.com/video_demo](https://www.evac-chair.com/video_demo)

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®. Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, EXCLUSIVE LIFETIME WARRANTY.

For further information please contact:

Evac+Chair North America LLC

3000 Marcus Ave, Suite #3E6, Lake Success, NY, 11042-1012

T 516-502-4240 W www.evac-chair.com E sales@evac-chair.com



THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR SINCE 1982