SUPER POWER EZ



Item# F-194PHBL-NA

The Super Power EZ Lift is an economical lift for above ground, backyard swimming pools and spas.

This lift is designed to clear walls up to 58 inches high. Push a button to raise or lower with 360 degrees of manual rotation.

The Super Power EZ Lift has all stainless steel construction and a durable powder coat finish.

+ FEATURES

- 400 lb (181 kg) weight apacity
- + Rechargeable battery with charger
- Clears walls up to 58" high (147 cm)
- Removable manual rotation bar
- Durable mesh sling- standard size: 38" x 20" (96 cm x 51 cm)
- Adjustable stainless steel chains
- 5 year structural warranty
- 5 year pro-rated electronics warranty

OPTIONS/ACCESSORIES:

F-612TLA	Anchor, standard in-ground, 3.25" x 8" deep
F-440HBC-B	Cover, standard lockable, blue
F-440HBC-G	Cover, standard lockable, gray
F-440HBC-T	Cover, standard lockable, tan
F-0320HSO-4C	Hard seat option
F-037EZLS	Sling, extra long - 50″x 20″ (127 cm x 50 cm)
F-036RSS	Sling, standard with rigid insert & seat belt
F-360YSLSB	Sling, youth with seat belt - 34" x 18" (86 cm x 46 cm)



PROUDLY SOLD BY:

THE SUPER POWER EZ LIFT







To confirm specific application go to www.aquacreek.com for our Pool Lift Selector Become a member of our new Dealer and Architect page. Access to CAD blocks, CSI specs and sales materials.



THE LEADERS IN RECREATION, FITNESS, AND ABILITY!