

)rS Red

Seat Color Black Vinyl

Body Colors Blue Re

## COMPANION

4 wheel GC440C

The Companion four-wheel scooter offers the most leg room and foot room in its class with a modern, rugged design. Features include easy disassembly; an adjustable LED headlight; ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion four-wheel is rugged, dependable and an outstanding value for a full-time, four-wheel scooter.

S


N



## GC440C

## Companion 4 wheel Specifications

Companion 4 w	meer spo
Model	GC440C
Medicare Code	Pending
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5.3 mph
Operating Range 1	18.9 miles (U1)
Ground Clearance Under Dec	8.4 miles (40AH) k 4.5"
Ground Clearance to Transaxl	
Turning Radius	56.5"
	40AH optional)
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
	ynamic R Series
Dimensions:	,
Length	47.5"
Width	23.5"
Height (ground to top	
of back of seat)	37" - 39"
Ground to Top of Seat	24" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	18" - 20"
Number of Seat Height Position	
Size of Increments	1"
Front Axle to Rear Axle	35.5"
Weight of Unit: (assemble	d) 194 lbs (U1) 210 lbs (40AH)
Front Half	74 lbs
Rear Half	48 lbs
Seat w/Arms	26 lbs
Batteries 21 lbs ea. (U1), 2	
	olid/Low Profile
Front	10.75"
Rear	10.75"
Rear Anti-Tippers	3"
Standard Seat 3:	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16" 18" X 16"
Width x Depth Color	Black Vinyl
Optional Seats:	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black Vinyl
Optional Seats 2:	Captain's Seat
Back Height - no headrest	19"
Back Height with headrest	29"
Width x Depth	18" x 18"
Color	Black Vinyl
Seat Options:	Zidek villyl
Pneumatic Rise Seat	Yes
Power Elevating Seat	Yes
Warranty:	. 63
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years
Batteries	1 Year





Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment



Patented infinite adjustable tiller



Quiet & maintenance free motor

## Standard Features on Every Model:

- Under mat troubleshooting guide
- Delta tiller with auto lock mechanism
- Easy reach tiller adjustment
- Adjustable sliding seat
- Adjustable LED headlight
- U-1 batteries standard





Disassembles easily!

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

