



FEATURES

- Simple push button remote
- Non AC Powered
- Whisper quite operation
- Fits standard bathtubs and 40" showers
- Adaptable for left or right hand tubs
- No plubming or structural changes needed
- Easilty portable to new locations
- Lifetime / 4 Year Electrical Warranty

SPECIFICATIONS

Lifting capacity:	350 lbs.
Seat Width:	15.75″
Seat Depth:	15.5″
Seat Height (Raised):	36″
Backrest Width:	15″
Backrest Height:	15″
Weight, Unit	95 lbs.
Reclines:	No
Power Source:	
18 Volt DC Craftsman Battery - UL Approved System	

DESCRIPTION

The USM Luxury Bath Lift is an exceptional product that will allow you to avoid having to step over the ledge of your tub. Unlike other lifts that require you to enter the tub and lower yourself to the chair height, the Luxury Bath Lift will raise you up, allowing for plenty of leg room to safely clear the edge of the tub, before safely placing you in your tub for a relaxing and therapeutic bath. This entire process takes just one minute and is easily repeated with the push of a button when you are ready to exit your tub. You can eliminate the risk of falling and safely enjoy your existing tub again with USM's Luxury Bath Lift.

Optiona Peri-Care Chair

With a bathtub lift from USM, you can easily care for a family member in your home or you can help give a family member an added security so they can maintain their independence.