

 \pm

S

 \simeq

 \subseteq

5



3 wheel GB107

This Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter, however, it has a weight capacity of 300 lbs. and increased foot room. The new Buzz LT also has a charger port in the delta tiller, a new hard plastic rugged basket, updated black wheels and a plush, padded seat. Additionally, the Buzz LT comes with a lifetime warranty on the frame, a two year warranty on the drive train and electronics, and comes standard with a one year warranty on the batteries. Enjoy its tight indoor maneuverability and smooth outdoor use with 3" of ground clearance!





3 wheel Specifications

*	
Model	GB107
Medicare Code	K0800 PENDING
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range ¹	8 miles
Ground Clearance Under Deck	3.0"
Turning Radius	37"
Type Batteries	2-12AH
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Dimensions:	
Length	40.5"
Width	21.5"
Height (ground to top of seat back)	33" - 35"
Ground to Top of Seat	21" - 23"
Ground to Top of Deck	4.75"
Deck to Top of Seat	16" - 18"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	31.5"
Weight of Unit: (assembled)	98 lbs
Front Half	33 lbs
Rear Half	25 lbs
Seat w/Arms	21 lbs
Batteries (pack with 2 batteries)	18 lbs
Tires:	Solid, Low Profile
Front	8"
Rear	8"
Rear Anti-Tippers	2"
Standard Seat:	Stadium Style
Back Height - no headrest	16"
Width x Depth	17" x 16"
Color	Black Vinyl
Warranty:	,
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Battery	1 Year
	i icui

- 1. Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.
- ^{2.} Electronics warranty excludes batteries.

Features

- Changeable color through chip resistant plastics
- Charger port in the tiller
- An increased weight capacity of 300 lbs.
- Increased foot room
- Chip resistant color through wheels





Quickly disassembles into five, lightweight pieces for easy transport!

