

# Roscoe Mobility Knee Scooter

The Knee Scooter offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Scooter is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required. The Knee Scooter is equipped with an easy folding mechanism that collapses using a thumb release lever and fits beneath the raised knee platform for ease of storage and transportation.

Black Finish



White Finish



NEW!  
Burgundy Finish



## Product Features

- Scooter can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Large 8-inch clear PU wheels with sealed bearings
- Adjustable hand brake ensures safety
- Padded knee platform with easy-adjustment mechanism
- Padded handlebar grips
- Limited lifetime warranty on frame, 90 days on all other parts

## Product Specifications

- Knee Platform Measurements: 7" (W) x 14" (L)
- Frame Finish: powder-coated black, white or burgundy
- Adjustable Handle Height: 31" - 41"
- Knee Pad Height Adjustment: 17 1/2" - 22"
- Product Weight: 22.8 lb.
- Weight Capacity: 350 lb.



**Knee Scooter  
Replacement  
Parts Available!**

### New and Improved Maintenance Free Axle and Axle Bearings!



Old Rear Axle

**NEW** Rear Axle



Adjustable Hand Brake



Handle Bar can be Adjusted

Item #	Description	Patient Height Range	Color	Qty
ROS-KSB	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	Black	each
ROS-KSW	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	White	each
ROS-KSBG	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	Burgundy	each



**To Order:**  
(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com