

Fingertip Pulse Oximeter

LED DISPLAY



The **Roscoe OTC Fingertip Pulse Oximeter** allows you to monitor and maintain pulse and oxygen saturation levels on-the-go. Accurate and fast oxygen level and heart rate measuring technology is ideal for sports enthusiasts during strenuous activities or mountain climbing. Recommended for spot-check monitoring only. No prescription required.



PRODUCT FEATURES:

- Approved for over-the-counter retail sale
- Large easy-to-read display
- Soft silicone padding for a comfortable fit
- One-button operation
- Up to 30 hours of use with batteries
- Auto-off function
- Low battery indicator
- Manual, lanyard and two AAA alkaline batteries included
- Attractive retail packaging
- One-year limited warranty

PRODUCT SPECIFICATIONS:

Dimensions:	2.25"L x 1.25"W x 1.25"H
Weight:	1.6 oz (including batteries)
Display Type:	LED
Display Range:	0 - 99%
Measurement Range:	0 - 100%
Accuracy:	70% - 100% ±3%
Resolution:	±2% for SpO2
Operating Temperature:	41-105°F (5~40°C)
Storage Temperature:	-4 -131°F (-20~55°C)
Ambient Humidity:	15%-80% in operation; ≤93% in storage

Item #	Description	UPC	UOM
POX-ROS	Roscoe Fingertip Pulse Oximeter	092237621981	each