



Infrared Forehead Thermometer

The Sejoy Infrared Thermometer delivers instant, single second temperature reads on both foreheads and objects alike. With a fever alarm and memory function to store up to 10 temperature readings, the Infrared Thermometer will be your go-to for quick and reliable temperature measurements. Dual scale readings in Fahrenheit and Celsius, low battery indicator and auto shut-off features also included. Retail package comes complete with Infrared Thermometer, 2 AAA Alkaline Batteries and User Manual.

FEATURES

- Fast 1 second readings
- Memory of up to 10 temperature readings
- Optional backlight for easy reading
- Results indicator lets you know when reading is ready
- Non-touch measurement for forehead or objects such as baby bottles
- Dual scale modes – Fahrenheit or Celsius
- Low battery indicator
- Fever alarm
- Auto shut-off

SPECIFICATIONS

| | |
|------------------------|--|
| Product Dimensions: | 6" L x 1.6" W x 1.7" D |
| Retail Box Dimensions: | 6.93" L x 3.66" W x 2.05" D |
| Product Weight: | 2.96 oz (with batteries) |
| Measure Time: | Approximately 1 second |
| Measuring Range: | Forehead Mode: 93.2°F ~ 109.4°F (34°C ~ 43°C) |
| | Object Mode: 32°F ~ 212°F (0°C ~ 100°C) |
| Battery Size: | AAA Battery (2) |
| Battery Life: | Approximately 1 year/6000 readings |
| Expected Service Life: | 3 Years |
| Warranty: | One year limited |

| Item # | Description | GTIN | UOM |
|---------|-------------------------------|----------------|------|
| DET-306 | Infrared Forehead Thermometer | 06970392215631 | each |