

PORTABLE PRO POOL 2 LIFT

US PATENT NO'S: [7,249,386 B2] [D507,769 S]




MEASUREMENT TAKEN WITH THE LIFT AT FULL LOAD



PRODUCT WEIGHT (NO COUNTERWEIGHT): 240lbs [108kg] PRODUCT WEIGHT (WITH COUNTERWEIGHT): 840lbs [381kg]

MAX WEIGHT CAPACITY: 350lbs [160kg]





9889 GARRYMORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602

NOTE:

*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

PRODUCT: PORTABLE PRO POOL 2		
PART #: F-PPP2		
SCALE:	1:24	DATE: 7/23/2020
DRAWN BY:	KR	REVISION: 01