


AL300-FULL

EXTERIOR HITCH MOUNT PLATFORM LIFT



A FLEXIBLE LIFT FOR SCOOTERS & POWER WHEELCHAIRS

The simple hitch-based installation requires little to no modification to your vehicle and allows for easy removal. Unmatched flexibility for users who expect to haul a different type of mobility device over time.

Consumers have the choice of a hold-down arm (power wheelchairs), hold-down foot (scooters), or retractable strap kit.



Hold-down arm
(Power wheelchairs only)



Hold-down foot
(Scooters only)



Retractable strap kit

AL300-FULL PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	- Power wheelchairs - Scooters
HITCH CLASS	II or greater
POWERED	Fully powered
AVAILABLE OPTIONS	- Battery Pack - Cover - License plate riser - Swing-away

INSTALLED WEIGHT	73 lb
PLATFORM SIZE	30.25" x 42"
WARRANTY	3-year
INCLUDED	- License plate mount - Wheel chocks
SECUREMENT OPTIONS	- Hold-down arm - Hold-down foot - Retractable strap kit



Harmar.com
800.833.0478

18MAR2022 | MKT-000187 A

