



# UNLEASH YOUR **ZOOM!**

NEW CAMO COLOR OPTIONS ARE HERE



CAPACITY



TOP SPEED



SUSPENSION



**BAJA** MOBILITY



[bajamobility.com](http://bajamobility.com)

## LIGHTER, FASTER, AND BUILT FOR THE BOLD.

Venture with confidence knowing you're equipped with the perfect blend of power and agility for all your journeys. An impressive ground clearance will get you through those rough spots with ease, and you'll be flying on flat terrain with a top speed of 8 mph. With a 400-lb. weight capacity and a lithium-ion battery pack that delivers up to a 29-mile range, the Bandit is built to defy the elements – and keep you moving all day.



# BAJA® BANDIT (BA140)



FDA Class II Medical Device\*

TAN  
(MATTE)

SAGE  
(MATTE)

CAMO

BLACK  
CAMO

WHITE  
CAMO



K0899 Power Mobility Device not coded by DME PDAC. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

\* The Baja® Bandit is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

1. Maximum speed varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of +10% or -5%.
2. Range per charge is established by testing conducted at maximum weight capacity, in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. These results occur from theoretical calculations based on battery specifications and drive system performance.
3. Turning radius can be subject to a variance of (+ or -) 3% due to continual product improvement and shifts in manufacturing tolerances.

## FEATURES

- Maximum speed of up to 8 mph
- 4.5" of ground clearance to tackle uneven terrain
- 400 lb. weight capacity
- Combine durability and speed and cruise up to 8 mph
- Competitive 53" turning radius
- Infinite tiller adjustment lever
- Convenient USB port for charging your devices
- Full LED lighting package for maximum visibility
- Full bumper package (standard)
- Supplemental hand brake for additional stopping power.
- Feather-touch disassembly with the heaviest piece (front) weighing 71 lbs.
- Seat post designed with full suspension for comfort and stability
- The upgradable dual lithium-ion battery configuration will deliver up to a 29-mile range on a single charge (14.5 miles with the standard single lithium-ion battery)
- 3-year warranty on lithium-ion battery

## ACCESSORIES

Only one rear accessory can be used at a time.

- 3" single or double cane/crutch holder
- 50/60/70" lap belt
- Cup holder
- Cell phone holder
- Forearm crutch holder
- Military patches
- Oxygen holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder

Cup holder, rear basket and saddlebag can be ordered as a bundle.

## SPECIFICATIONS

WEIGHT CAPACITY	400 lbs.	
MAXIMUM SPEED <sup>1</sup>	Up to 8 mph	
DIMENSIONS <sup>3</sup>	Width: 27.83"	Length: 49.6"
TURNING RADIUS <sup>3</sup>	53"	
RANGE PER CHARGE <sup>1,2</sup>	<b>20Ah Battery Pack:</b> Up to 14.5 miles at 220 lbs. Up to 11.25 miles at 400 lbs.  <b>40Ah Battery Pack:</b> Up to 29 miles at 220 lbs. Up to 22.5 miles at 400 lbs.	
TOTAL WEIGHT	183.4 lbs. with 20 Ah battery pack 197.8 lbs. with 40 Ah battery pack	
HEAVIEST PIECE	71 lbs. (front)	

## WARRANTY

FRAME: Lifetime limited  
ELECTRONICS: 2-year limited

DRIVETRAIN: 2-year limited  
BATTERY: 3 years

AUTHORIZED PROVIDER

bajamobility.com

