

Rehab Line Heavy Duty Beds



Bed Features	Night Rider HD	Night Rider SHD	Valiant HD	Valiant SHD	Companion HD	Companion SHD
Deck Height	13.5"	13.5"	13.5"	13.5"	13.5"	13.5"
Range	10"	10"	10"	10"	N/A	N/A
Head & Foot	✓	✓	✓	✓	✓	✓
Hi-Low	✓	✓	✓	✓	✗	✗
Trendelenburg	✗	✗	✓	✓	✗	✗
Reverse Trendelenburg	✗	✗	✓	✓	✗	✗
Cardiac Chair	✗	✗	✓	✓	✗	✗
Twin 80 & 84	✓	✓	✓	✓	✓	✓
Full 80 & 84	✓	✓	✓	✓	✓	✓
Dual King	✓	✓	✓	✓	✓	✓
Special Size 48x80	✓	✓	✓	✓	✓	✓
Capacity	600 lbs	750 lbs	600 lbs	750 lbs	600 lbs	750 lbs
Standard Hand Control	Wired	Wired	Wired	Wired	Wired	Wired
Deluxe Upgrade Available	✓	✓	✓	✓	✓	✓
Hoyer Compatible*	Low Profile	Low Profile	Low Profile	Low Profile	Low Profile	Low Profile
Environmental Control Unit (ECU) Support	✓	✓	✓	✓	✓	✓

Technical Specs	Night Rider HD	Night Rider SHD	Valiant HD	Valiant SHD	Companion HD	Companion SHD
Head Angle	60°	60°	60°	60°	60°	60°
Foot Angle	35°	35°	35°	35°	35°	35°
Trendelenburg Angle	N/A	N/A	10°	10°	N/A	N/A
Reverse Trendelenburg Angle	N/A	N/A	10°	10°	N/A	N/A
Cardiac Chair Angle	N/A	N/A	70°	70°	N/A	N/A
Twin 80 & 84 Frame Weight**	275 lbs	285 lbs	275 lbs	285 lbs	227 lbs	237 lbs
Full 80 & 84 Frame Weight**	323 lbs	333 lbs	323 lbs	333 lbs	266 lbs	276 lbs
48x80 Frame Weight**	298 lbs	308 lbs	298 lbs	308 lbs	266 lbs	276 lbs
Queen Frame Weight**	332 lbs	342 lbs	332 lbs	342 lbs	275 lbs	285 lbs
Current	AC	AC	AC	AC	AC	AC
Duty Cycle	20%	20%	20%	AC	AC	AC

** Crated, without mattress